




Ambu® Neuroline Monopolar

Single-patient EMG Needle Electrode



The Ambu Neuroline Monopolar promises to deliver excellent control and flexibility with every use.

SUITABLE APPLICATION

Electromyography (EMG)

FEATURES

- The tip is superior in its sharpness for low penetration resistance
- The needles have an ultra-smooth coating, which reduces friction when inserting the needle and

significantly minimises patient discomfort when working in the target area

- Made of stainless steel and carefully manufactured to ensure the optimum combination of strength and flexibility
- The ribbed and ergonomically shaped hub design allows for easy and secure handling
- Choice of five convenient sizes - all colour-coded
- Available with or without a 76 cm (30") pre-attached lead wire

ORDERING INFORMATION

NEEDLES WITHOUT LEAD WIRES

| ITEM NO. | LENGTH | CALIBER | HUB COLOUR | CONNECTOR |
|--------------|--------|------------|---------------|------------------|
| 742 25-36/40 | 25 mm | 0.36 (28G) | ■ Grey | 1.5 (DIN 42 802) |
| 742 38-36/40 | 38 mm | 0.36 (28G) | ■ Light green | 1.5 (DIN 42 802) |
| 742 38-45/40 | 38 mm | 0.45 (26G) | ■ Green | 1.5 (DIN 42 802) |
| 742 50-45/40 | 50 mm | 0.45 (26G) | ■ Yellow | 1.5 (DIN 42 802) |
| 742 75-50/40 | 75 mm | 0.50 (25G) | ■ Red | 1.5 (DIN 42 802) |

NEEDLES WITH LEAD WIRES

| ITEM NO. | LENGTH | CALIBER | HUB COLOUR | CONNECTOR |
|--------------|--------|------------|---------------|------------------|
| 743 25-36/40 | 25 mm | 0.36 (28G) | ■ Grey | 1.5 (DIN 42 802) |
| 743 38-36/40 | 38 mm | 0.36 (28G) | ■ Light green | 1.5 (DIN 42 802) |
| 743 38-45/40 | 38 mm | 0.45 (26G) | ■ Green | 1.5 (DIN 42 802) |
| 743 50-45/40 | 50 mm | 0.45 (26G) | ■ Yellow | 1.5 (DIN 42 802) |
| 743 75-50/40 | 75 mm | 0.50 (25G) | ■ Red | 1.5 (DIN 42 802) |

Please see page 22 for cables