

# BIOCRYO

COLD COMPRESSION  
THERAPY SYSTEM

**Recover Faster**  
Cold Compression Made Easy



GC-3045 BioCryo Full Leg  
(Available in S/M/L)

# BIOCRYO


## COLD COMPRESSION THERAPY SYSTEM

### How Does the BioCryo System Work?

1. The BioCryo starts with sequential compression that massages the tissue to move fluids through and out of the injured area or surgical site, speeding the healing process and preventing infection.



2. The BioCryo combines this massaging action with the convenient BioCryo gel pack insert which conforms to the body for maximum surface contact, pulling heat away from the injury site while reducing inflammation and pain.



### What are the BioCryo advantages?

- Convenient cold gel pack insert eliminates the need for ice and water.
- 4 sequential chambers of cold compression therapy is more effective than single chamber devices.
- Pain management without the use of narcotics helps patients reduce their need for prescription medication.
- Increased patient compliance- easy to use and comfortable to wear.
- Shorten post surgical recovery time.
- Light weight and portable, the BioCryo is easy to carry and set up.
- Gel packs store conveniently in the freezer and stay ice cold for hours of use.
- Unilateral and/or Bilateral functionality.
- Wide selection of sleeves including latex free custom garments for unique sizing.
- Medically proven results to prevent infection and relieve swelling, soreness, and fatigue.

The BioCryo gel  
insert easily  
attaches to all  
compression  
sleeves





# The BioCryo Gel Insert

Our BioCryo gel inserts are latex free and made of high quality 200 Denier coated oxford nylon fabric with 3 mils of polyurethane. These uniquely designed inserts flex and conform to any extremity, delivering cold compression without having to add or pump ice water to the garment. These inserts are to be used inside our standard, adjustable, and custom garments.



BioCryo Gel Insert CP-S/M/L



BioCryo Shoulder GC-3035-SH-C



BioCryo Back GC-3035-B



BioCryo Arm & Shoulder GC-3035-SH



BioCryo Elbow GC-3035-EC



BioCryo Ankle GC-3045-A



BioCryo Full Arm GC-3035-L (Available in S/M/L)



BioCryo Knee GC-3045-KC



BioCryo Hand/Wrist GC-3035-H-W

# BIOCRYO

## Cold Compression Therapy from the #1 Compression Device Maker

Our cold compression therapy system combines active cold and sequential compression therapies to provide maximum power to manage pain, reduce inflammation, and speed healing for a wide variety of ailments. BioCryo technology stimulates tissue repair by enhancing lymphatic function, removing edema from the injury site, and encouraging oxygenated blood flow to the site, leading to a faster recovery. Get back in action with BioCryo from Bio Compression Systems, the leading manufacturer of Pneumatic Compression Devices.



SC-2004-C

### Contraindications:

Use of this device is contraindicated for patients with the following conditions:

- Congestive heart failure • Deep Vein Thrombosis
- Inflammatory Phlebitis or episodes of Pulmonary Embolism
- Infections in the limb, including cellulitis without appropriate antibiotic coverage
- Presence of cancer unless for palliative care • Cold Urticaria • Cold Erythema
- Cold Hemoglobinuria • Other indications as identified by the treating physician



Bio Compression Systems  
120 West Commercial Avenue  
Moonachie, NJ 07074  
Phone: 800-888-0908/ Fax: 201-939-4503  
email: biosystems@biocompression.com  
www.biocompression.com