


# Imaging and Smart Dilation™ solutions for the alimentary canal


endoflip®  
Endolumenal Functional Lumen Imaging Probe



## Esophageal Smart Dilation™ Catheter

DESCRIPTION	PART NUMBER
endoflip® System	EF-100
endoflip® Cart	DD-360
endoflip® 8cm measurement catheter with pressure	EF-325N
endoflip® 16cm measurement catheter with pressure	EF-322N
endoflip® 20mm esophageal Smart Dilation™ catheter	ES-320
endoflip® 30mm achalasia Smart Dilation™ catheter	ES-330

# Patient with dysphagia/food impaction



- Smart dilator – measurement nodes display stricture size before and after balloon dilation
- Better control – hands-free dilation with the endoflip® system
- No fluoroscopy required – removing exposure of patient and staff to radiation during dilation
- Precision – only one balloon required to read esophageal diameter directly from the endoflip® screen