

EMG ELECTRODES

Ground electrodes and connecting cables



32ERD81

Ground plate electrode

diameter 31mm, 1.5mm DIN
connector, 1200mm long

1 PU = 1 pc

Price list page 3



32KABEL9

**Snap connector cable for
ground electrodes**

with 4mm connector (or with
2mm connector) suitable for
ground electrodes 32ERD5 ...
32ERD71, 1500mm long

1 PU = 1 pc

Price list page 5



32KABEL9A

**Snap connector cable for
ground electrodes**

with DIN connector, suitable for
ground electrodes 32ERD5 ...
32ERD71, 1500mm long

1 PU = 1 pc

Price list page 5



32KABEL0010

**Snap connector cable for
ground electrodes**

with DIN connector, suitable for
32ERD0010, 32ERD0011 and
32ERD0012, 1,500mm long

1 PU = 1 pc

Price list page 4

EMG ELECTRODES

Ground tape



32ERD5

Silicone ground tape for children

12mm x 200mm (to be used in dry condition), colour: black

1 PU = 1 pc
Price list page 3



32ERD51

Copper ground tape for children

20mm x 250mm (to be moistened before use), colour: green

1 PU = 1 pc
Price list page 3



32ERD52

Copper ground tape

20mm x 350mm (to be moistened before use), colour: green

1 PU = 1 pc
Price list page 3



32ERD6

Silicone ground tape for adults

20mm x 450mm (to be used in dry condition), colour: black

1 PU = 1 pc
Price list page 3



32ERD61

Copper ground tape for adults

20mm x 450mm (to be moistened before use), colour: green

1 PU = 1 pc
Price list page 3



32ERD7

Silicone ground tape for adults

20mm x 700mm (to be used in dry condition), colour: black

1 PU = 1 pc
Price list page 3



32ERD71

Copper ground tape for adults

20mm x 700mm (to be moistened before use),
colour: green

1 PU = 1 pc
Price list page 3



32ERD8

Copper ground tape for adults

with preassembled 1500mm
cable, DIN connector, 20mm x
550mm long ground strap with
fixing tape

1 PU = 1 pc
Price list page 3



32ERD0010

Elastic ground tape

with snap connector,
400mm long, colour: blue

1 PU = 1 pc
Price list page 3



32ERD0011

Elastic ground tape

with snap connector,
250mm long, colour: blue

1 PU = 1 pc
Price list page 3



32ERD0012

Elastic ground tape

with snap connector,
700mm long, colour: blue

1 PU = 1 pc
Price list page 3

EMG ELECTRODES

Adhesive electrodes with connection cables



32KE7

**Ag/AgCl adhesive electrode
GVB**

diameter 10mm, 1000mm
connection cable with DIN
connector

1 PU = 1 pc
Price list page 6



32KE13

**Disposable adhesive electrodes
"SENSITABS"**

20mm x 20mm contact area,
without cables

1 PU = 100 pcs
Price list page 5



32KE13-1

**Disposable adhesive electrodes
"MAX-TAB"**

25mm x 13mm contact area,
without cables

1 PU = 100 pcs
Price list page 5



32KABEL17

Pair of alligator clips

for disposable adhesive
electrodes (1 pair), with 1.5mm
DIN connector, matching
connection cable for 32KE13,
1000mm long

1 PU = 2 pcs
Price list page 4



32KA0013

Pair of alligator clips

for disposable adhesive
electrodes, for connection to
multi-plug system, 150mm long

1 PU = 1 pc
Price list page 4

EMG ELECTRODES

Needle connection cable



32KA10

Needle connection cable NK

for reusable EMG needle electrodes NK, 1000mm long

1 PU = 1 pc
Price list page 4



32KA20

Needle connection cable SCHULER

for disposable EMG needle electrodes, suitable for 32NE69-1 to 32NE73-1 and NEUROWERK EMG devices as of 08/2015 (serial number 5690), 1200mm long

1 PU = 1 pc
Price list page 4



32KA20A

Needle connection cable SCHULER

for disposable EMG needle electrodes, suitable for 32NE69-1 to 32NE73-1 and NEUROWERK EMG devices up to 07/2015, 1200mm long

1 PU = 1 pc
Price list page 4



32KA0010

Needle connection cable

for "Myoline" needle electrodes, multi-plug system, 5-pin 270° connector, suitable for NEUROWERK EMG devices up to 07/2015, 1200mm long

1 PU = 1 pc
Price list page 4





32KABEL95-2

Connection cable for monopolar needle electrodes

for connection to 1-channel electrode cable, 1.25mm female connector to DIN connector, 500mm long

1 PU = 3 pcs

Price list page 5



32KABEL95-3

Connection cable for monopolar needle electrodes

for connection to 1-channel electrode cable, 1.25mm female connector to DIN connector, 1000mm long

1 PU = 3 pcs

Price list page 5



EMG ELECTRODES

Surface electrodes



32AE4

Surface electrode NK

5 electrodes and 1 jumper cable

1 PU = 5+1 pcs

Price list page 2



32AE4-1

Jumper cable NK – single

1 PU = 1pc

Price list page 2




32AE5

Surface electrode NK

sintered Ag/AgCl surface
electrode

1 PU = 1 pc

Price list page 2



32AE6-1

Jumper cable NK – single

with 2 x 1.5mm DIN connector

1 PU = 1 pc

Price list page 2



32AE7

Surface electrode SCHULER

240° 5-pin DIN circular connector, 1200mm long

1 PU = 1 pc
Price list page 2



32AE7-1

Surface electrode SCHULER


with DIN connector, 1 pair, suitable for 32KABEL7, 200mm long

1 PU = 2 pcs
Price list page 2



EMG ELECTRODES

Surface stimulation electrodes



32AE1-S

Block electrode

5-pin DIN connector, incl. 10 felts and 1 fixing tape

1 PU = 1 pc
Price list page 2



32RE1

Surface stimulation electrode NK for children

incl. 4 round felts and fixing tape, connection cable 2000mm long, circular connector 240°

1 PU = 1 pc
Price list page 9



32RE2

Surface stimulation electrode NK for adults

incl. 4 round felts and fixing tape, connection cable 2000mm long, circular connector 240°

1 PU = 1 pc
Price list page 10



32RE3

Surface stimulation electrode GVB for adults

incl. 4 round felts and fixing tape, connection cable 1500mm long, circular connector 240°

1 PU = 1 pc
Price list page 10



32RE7

**Surface stimulation electrode
SCHULER for adults**

3-pin DIN connector 180°,
cable 2000mm long, suitable for
"NeuroScreen"

1 PU = 1 pc
Price list page 10



32RE7-1

**Surface stimulation electrode
SCHULER for adults**

5-pin DIN connector 240°,
cable 2000mm long

1 PU = 1 pc
Price list page 10



32RE7-2

**Surface stimulation electrode
SCHULER for adults**

2 x 1.5mm DIN connectors,
cable 2000mm long

1 PU = 1 pc
Price list page 10



32RE4

Loop electrode NK

with preassembled cable,
2000mm long, circular
connector 240°

1 PU = 1 pc
Price list page 10





32RE6

Loop electrode SCHULER

with 1 replacement loop;
loops interchangeable separately,
cable 1200mm long, circular
connector 240°

1 PU = 1 pc
Price list page 10



32RE6-1

**Replacement loops for loop
electrode 32RE6**

diameter 2.6mm, pluggable

1 PU = 2 pcs
Price list page 10



32RE6-11

**Replacement loops for loop
electrode 32RE6**

diameter 0.7mm, pluggable

1 PU = 2 pcs
Price list page 10



32RE6-3

Replacement connection cable

DIN connector, suitable for
32RE6, 1200mm long

1 PU = 1 pc
Price list page 10



32RE7-3

Wire loop electrode GVB (thin)

suitable for multi-function
cables, 150mm long

1 PU = 1 pc

Price list page 10



32RE0010

**Finger-toe stimulation
electrode GVB**

suitable for multi-function
cables, 150mm long

1 PU = 1 pc

Price list page 9



32KA0011

Multi-function cable GVB

connection cable for multi-plug
system (mot. and sens. neurography),
5-pin DIN connector 240°, suitable
for NEUROWERK EMG devices up to
07/2015, 1200mm long

1 PU = 1 pc

Price list page 4



32KA0012

Multi-function cable GVB


connection cable for multi-plug
system, 5-pin DIN connector
240°, suitable for NEUROWERK
EMG devices as of 08/2015 (serial
number 5690), 1200mm long

1 PU = 1 pc

Price list page 4

EMG ELECTRODES

Accessories for surface stimulation electrodes



32MSE-F

**Reusable metal tips SCHULER
for 32RE7 and 32RE7-1**

gold-plated brass, flat design

1 PU = 2 pcs
Price list page 8



32MSE-R

**Reusable metal tips SCHULER
for 32RE7 and 32RE7-1**

gold-plated brass, round design

1 PU = 2 pcs
Price list page 8



32RE2-M

Reusable metal tips NK
suitable for stimulation electrode
32RE2

1 PU = 1 pc
Price list page 10



32BAND4

Fixing tape with hole
suitable for surface stimulation
electrodes 32RE7 and 32RE7-1,
420mm long

1 PU = 1 pc
Price list page 2



32BAND6

Fixing tape
suitable for block electrode
32AE1-S, 450mm long

1 PU = 1 pc
Price list page 2



32FILZ1-S

Rectangular replacement felts

suitable for block electrode
32AE1-S

1 PU = 10 pcs
Price list page 3



32FILZ2

Round replacement felts

suitable for surface stimulation
electrode 32RE1, diameter 3mm,
8mm long

1 PU = 10 pcs
Price list page 3



32FILZ3

Round replacement felts

suitable for surface stimulation
electrodes 32RE2 and 32RE3,
diameter 7mm, 8mm long

1 PU = 10 pcs
Price list page 3



32FILZ4

Round replacement felts

suitable for surface stimulation
electrodes 32RE7 and 32RE7-1,
diameter 6mm, 7mm long

1 PU = 20 pcs
Price list page 3



EMG ACCESSORIES

Other EMG accessories



32KABEL7

1-channel electrode cable

with 5-pin DIN connector 240°
to 3 x 1.5mm DIN socket,
1500mm long

1 PU = 1 pc
Price list page 4



36HTS409B

Skin temperature sensor

to measure the skin surface
temperature (only suitable for
NEUROWERK EMG), 3000mm
long

1 PU = 1 pc
Price list page 12